


Sample & Assay Technologies **Microbiome: From Identification to Characterization**




Ivorlyne Greene Ph.D.
QIAGEN Inc.

Sample & Assay Technologies **Legal Disclaimer**

- QIAGEN products shown here are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.
- For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.QIAGEN.com or can be requested from QIAGEN Technical Services or your local distributor.

Title Location Date 2

Sample & Assay Technologies **Agenda**



- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Title Location Date 3

Sample & Assay Technologies

Agenda


- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Titlu, Location, Date

Sample & Assay Technologies

Humans or Superorganisms?

Cellular composition of the organism



Human: Estimations of the number of microbial cells that live in and on the human body, human cells are outnumbered by a factor of 10.

Nomenclature:

- **Microbiota** are the microbes that live in a specific location, e.g. the human body, the gut, soil, etc.
- **Metagenomics** is the study of the collection of genomes derived from a specific sample or community
- **Microbes** are microscopic organisms that can be either single or multicellular.


For the purposes of this webinar, we will use **microbiome** to describe the collective genomes of the microbiota that inhabit a specific location.

Titlu, Location, Date

Sample & Assay Technologies

Microbiota composition

Microorganisms cluster by body site



Cataloguing efforts by the NIH Human Microbiome Project suggest:

- ~10,000 organisms live with us
- ~ 8×10^6 genes in this "second genome"

Identifying microbiota in healthy individuals revealed:

- Different body sites have unique communities
- Race, Age, Gender, Weight or Ethnicity have an effect

Titlu, Location, Date

Sample & Assay Technologies **Complexity and function of genomic content**

Function of microbiome enables individual survival

Each organism has developed genetic content for its own survival in a specific environment

Metabolism tuned to local nutrient sources

Virulence factors for stable colonization

Antibiotic resistance genes to metabolize toxins

81231_521169_High-Resolution Microarray

Title, Location, Date 7

Sample & Assay Technologies **Agenda**

- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Title, Location, Date 8

Sample & Assay Technologies **Current Methods for microbial analysis**

Identifying microbes and their genes

- Culture
 - Nutrients
 - Antibiotics
- 16S rRNA gene clone library construction
 - Pan 16S rRNA PCR amplification → cloned → sequenced
- Microarray
- Next generation sequencing
 - 16S rRNA sequencing
 - Whole genome sequencing
- MALDI
- qPCR
 - Target dependent (16S rRNA or gene)

Title, Location, Date 9

Sample & Assay Technologies **16S rRNA gene as a phylogenetic marker for bacterial ID**

16S rRNA gene

1 542

■ - Conserved region ■ - Variable region

- Classification from the variable sequences
- 16s rRNA sequence similarity

95% genus level, **97% species level**, 99% strain level

Titb, Location, Date 10

Sample & Assay Technologies **How did we get started?**

First report of the composition of human body site

Science VOL 312 2 JUNE 2006

Metagenomic Analysis of the Human Distal Gut Microbiome

Steven R. Gill,^{1,2*} Mihai Pop,^{1†} Robert T. DeBoy,¹ Paul B. Eckburg,^{2,3,4} Peter J. Turnbaugh,⁵ Buck S. Samuel,⁵ Jeffrey I. Gordon,⁵ David A. Relman,^{2,3,4} Claire M. Fraser-Liggett,^{3,4} Karen E. Nelson*

If we can do one....can we do them all....

Titb, Location, Date 11

Sample & Assay Technologies **International efforts to catalog the microbiome**


HMP Meta-HIT

International Human Microbiome Consortium

Titb, Location, Date 12

Sample & Assay Technologies
Catalog by five primary body sites

NIH Human Microbiome Project



Profiling 5 body sites


- Nasal
- Mouth
- Skin
- Gastrointestinal system
- Urogenital

Compare between individuals:

- Healthy vs. Disease
- Treated vs. Untreated
- Twin studies
- Diet
- ...

Title, Location, Date 13

Sample & Assay Technologies
Agenda




- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Title, Location, Date 14

Sample & Assay Technologies
Association of the human microbiota and disease



- Gut
 - Intestinal infections
 - Obesity
 - Inflammatory Bowel Disease
- Airway
 - Pneumonia and other respiratory infections
 - Chronic Obstructive Pulmonary Disease
 - Cystic Fibrosis
- Urogenital
 - Bacterial Vaginosis
 - Urinary Tract Infections
 - Sexually Transmitted Disease
- Blood
 - Sepsis/Blood-stream infections
- Cancer
- Heart disease
- Neurological disorders
- Oral
 - Periodontitis
 - Gingivitis



Title, Location, Date 15

Sample & Assay Technologies **Physiological associations lead to new funding**

NIH funding across institutes for microbiome-related studies

2008

2013

Title, Location, Date 16

Sample & Assay Technologies **Focused Metagenomics Applications**

Examples from the next wave of microbiome experiments

- Screening for microbial genes in a specific metagenomic sample
Experiment 1: Antibiotic resistance genes in human gut
- Relationships between microorganisms that permit colonization
Experiment 2: Profile changes in vaginal flora during bacterial vaginosis
- Surveillance of samples for specific targets
Experiments 3: Screen samples for specific pathogens, resistance genes or virulence factors

Title, Location, Date 17





Sample & Assay Technologies **Antibiotic resistance gene reservoirs in the body**

Screening of the gut for presence of antibiotic genes



Gene	Resistance classification	1	2	3	4	5
aacC2	Aminoglycoside resistance					+
ermB	Macrolide lincosamide streptogramin B	+	+	+	+	+
mefA	Macrolide lincosamide streptogramin B	+	+	+	+	+
tetA	Tetracycline efflux pump	+	+	+		

500 ng of genomic DNA from stool samples originating from five healthy adults were analyzed for presence of antibiotic resistance genes. 87 unique antibiotic resistance genes were tested for by qPCR. Positive (+)/negative (blank) result for each antibiotic resistance gene was determined using a threshold cycle cutoff.

Title, Location, Date 18



Antibiotic resistance genes in our food supply?

Genes	Resistance classification	Sevage
aacC1	Aminoglycoside resistance	+
aadA1	Aminoglycoside resistance	+
GES	Class A beta-lactamase	+
SHV	Class A beta-lactamase	+/-
SHV (156G)	Class A beta-lactamase	+/-
SHV (238G/240E)	Class A beta-lactamase	+/-
T1A.1	Class A beta-lactamase	+
VEB	Class A beta-lactamase	+
ACT-1 group	Class C beta-lactamase	+/-
LAT	Class C beta-lactamase	+/-
MIR	Class C beta-lactamase	+/-
MCOX	Class C beta-lactamase	+
OXA-10 group	Class D beta-lactamase	+
OXA-2 group	Class D beta-lactamase	+
AAC(6)-Ib-cr	Aminoglycoside resistance	+
QnrB-5 group	Fluoroquinolone resistance	+
Qnr5	Fluoroquinolone resistance	+
ermB	Macrolide lincosamide streptogramin B	+
mefA	Macrolide lincosamide streptogramin B	+
tetA	Tetracycline efflux pump	+

Title, Location, Date 22

Common themes in these experiments


How best to carry out these focused metagenomic studies?

- Each case has the user profiling a distinct set of genes or species
- Experiments easily translate to large multi-factor experiments
- Range of sample types
- Quantitative and Qualitative data generation

Title, Location, Date 23

Agenda

- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions



Title, Location, Date 24

Sample & Assay Technologies
16S rRNA gene as a phylogenetic marker for bacterial ID

Sequencing or Real-time PCR (qPCR)

16S rRNA gene

■ - Conserved region ■ - Variable region

- Classification from the variable sequences
 - 16s rRNA sequence similarity
 - 95% genus level, **97% species level**, 99% strain level
- Assay design approach
 - use only sequences with taxonomy classified by the GreenGenes taxonomy
 - fairly specific probe + fairly specific primer pair = specific assay

Titb, Location, Date 25

Sample & Assay Technologies
How do hydrolysis probe assays work?

Exonuclease activity of Taq polymerase disrupts FRET by separating reporter from quencher

Titb, Location, Date 26


Sample & Assay Technologies
Performance testing of each assay

Dilution series testing for PCR efficiency and sensitivity


A Fluorescence (493-533)


B $y = -3.3236x + 40.236$
 $R^2 = 0.9965$


Titb, Location, Date 27




- Sample & Assay Technologies **QIAGEN's Microbial DNA qPCR assay pipeline**
- Develop an assay pipeline to support microbiome research
- Over 500 assays that target species-specific or gene-specific microbial DNA
 - >300 Bacteria identification assays
 - 8 Fungi identification assays
 - 1 Protist identification assay
 - 87 Antibiotic resistance genes
 - 87 Virulence factor genes
- Titlu, Locațiun, Data 29







Sample & Assay Technologies **Agenda**





- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions


Titlu, Locațiun, Data 33

Sample & Assay Technologies **Overview of Experimental Protocol**

Focused metagenomic experiments

Test for one gene or organism 
 Microbial DNA qPCR Assays
 Microbial DNA qPCR Assay Kits


Test for 2-6 genes or species 
 Microbial DNA qPCR Multi-Assay Kits

Test for 24-96 genes or species 
 Microbial DNA qPCR Arrays
 Custom Microbial DNA qPCR Arrays

Title, Location, Date 34

Sample & Assay Technologies **Isolation of microbial DNA from metagenomic samples**


Specialized kits and protocols to match samples





Sample Material	DNA Purification Kit	Catalog Number
Stool	QIAamp DNA Stool Mini Kit	51504
Blood	QIAamp UCP PurePathogen Blood Kit	50112
Cervical swab in transport media	QIAamp MinElute Media Kit	57414
Blood culture, bronchoalveolar lavage, carious dentine, cervical swab, isolated bacterial colony, sputum, saliva, swabs	QIAamp UCP Pathogen Mini Kit	50214


Title, Location, Date 35

Sample & Assay Technologies **Assay User Protocol**



Set up 4 PCR reactions per sample that include positive and negative controls. 

Run real-time PCR to obtain raw Ct values for each assay. 

Use analysis software to identify which species/gene are present. 

Title, Location, Date 36

Microbial DNA qPCR Array

Sample & Assay Technologies

Pre-printed assays profile up to 90 different species/genes

DNA Isolation → Detection by qPCR → Data Analysis

PCR plates (either 96-well or 384-well) are pre-printed with primers and probes.

Each numbered well is a separate assay that tests the same sample.

Integrated control assays:
 Host assays detect genomic DNA to test sample collection
 Pan A/C is a pan-fungal assay that detects the presence of fungal rRNA
 PanB1 and PanB2 detect bacterial 16S rRNA to determine bacterial load in the sample
 PPC is a positive PCR control reaction that tests if the PCR reactions failed from PCR inhibitors from the sample, etc.

Titlu, Location, Date

Layout of a Microbial DNA qPCR Array

Sample & Assay Technologies

Different arrays have different number of assays and samples

Titlu, Location, Date

Microbial DNA qPCR Arrays

Sample & Assay Technologies

Profile or identify the presence of microbial DNA with any array

Identification experiment answers the following question:
 Are any of these microbes or genes present in the sample?
 Must be compared against a known negative sample.
 Answers are Yes or No.

Profiling experiment answers the following question:
 Have the amounts of any of these microbes or genes changed?
 Must be compared against a reference sample.
 Answers are fold change.

Titlu, Location, Date

Sample & Assay Technologies **Microbial DNA qPCR Array Content**

14 Arrays

- Antibiotic Resistance Genes
- Bacterial Vaginosis
- Biodefense
- Food testing: Dairy
- Food testing: Meat
- Food testing: Poultry
- Food testing: Seafood
- Food testing: Vegetable
- Intestinal infections
- Respiratory Infections
- Sepsis
- Urinary Tract Infections
- Vaginal Flora
- Water Analysis

Titix Location Date 40

Sample & Assay Technologies **Microbial DNA qPCR Arrays Protocol**

Identification or Profiling

Isolate microbial genomic DNA using QiaAMP kit dependent upon sample.

Mix sample DNA with Microbial qPCR mastermix and aliquot into each well of plate. If performing for identification analysis, then run No Template Control sample.

Run real-time PCR to obtain raw Ct values for each assay.

Use analysis software to generate species profile or species/gene identification.

Profiling **Identification**

Titix Location Date 41

Sample & Assay Technologies **Application examples performed using these arrays**

Focused metagenomic studies to identify and profile

Vaginal Flora qPCR Array

Antibiotic Resistance Gene qPCR Array

Gene	Resistance classification	1	2	3	4	5
aacC2	Aminoglycoside resistance					+
ermB	Macrolide lincosamide streptogramin B	+	+	+	+	+
mefA	Macrolide lincosamide streptogramin B	+	+	+	+	+
tetA	Tetracycline efflux pump	+	+	+		

Titix Location Date 42

Sample & Assay Technologies

Agenda

- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Titix Location Date 43

Sample & Assay Technologies

Summary




- Defined the microbiome and the key role it plays in human health
- Discussed the focused metagenomic studies
- Introduced the Microbial DNA qPCR assay pipeline
 - Largest collection of Microbial ID assays available
 - Available in various product formats to simplify experiments

Titix Location Date 44

Sample & Assay Technologies

Agenda



- Humans or Super-organisms?
 - Introduction to the microbiome
- Cataloging our "Second Genome"
 - Technologies that unlocked metagenomics
- Microbiota: Biomarkers or Effectors?
 - Population Studies and disease associations
- Identify and profile relevant targets
 - How to design assays for the microbiome
- Integrating focused metagenomics into my research
 - Solutions for identification of profiling
- Summary
 - Identification to Characterization
- Questions

Titix Location Date 45

QIAGEN Sample & Assay Technologies


Questions?

Ask now or contact Technical Support

North America
Telephone: 888-503-3187
Email: support@qiagen.com

Europe
00800-22448000
Email: sab.eu@qiagen.com

Worldwide
Email: support@qiagen.com



Title Location Date
